


**MAS SD Recording System
Video + (Optional Sound)**



INTERFACE CABLES		
<i>Type</i>	<i>Connects</i>	<i>Purpose</i>
Cable A - Bespoke	Trigger to Control	System: Control of Recording Start / Pause
Cable B - Bespoke	Control to Recorder - REMOTE 1	System: Control of Recording Start / Pause
Cable C - Video / Composite BNC	X-ray to Isolator Isolator to Scan Converter Isolator to Recorder Recorder to Monitor	Video
Cable E - Audio / Phono to Jack	Pre Amp to Recorder Recorder to Monitor	Audio
Cable F - USB	Keyboard to Recorder - USB 3	Data