



Medical Limited

Dental / Surgical Loupes

Choice of Dental or Surgical Loupes involves Multiple Options to Consider

Comfortable Working Distance

Establish an optimal working distance from the patient.
Measure the distance between your Eyes and the Patient

Depth of Field

Viewing Distance Range that Loupes remain in focus
Provides ability to move closer / further away from the patient, without loss of Loupe focus

Dental Loupe Magnification

Options Range - 2.5x Magnification to 6.0x Magnification for most standard Surgical and Dental Loupe

We advise an initial Loupe with low magnification - 2.5x / 3.0x

Adjustment of Loupes is simpler with lower magnification

Field of View

Lower Magnification comes accords with larger Field of View

Higher Magnification reduces Field of View - recommended for those experienced in the use of Loupes, and confident with reduced Field of View

Clip-on Loupes

Clip-on Loupes are the most cost effective choice for those who already wear a compatible frame for attachment

Clip-on Loupes offer the same level of Magnification Quality as traditional Galilean Style Loupes, but with added flexibility

Clip-on Loupes Flip similar to Galilean Style, providing a normal view outwith the Loupe

Loupes are supplied with adjustment for the Pupillary Distance of Loupe

Clip-on Loupes can be the most flexible and effective solution for Surgeons, Dentists, Dental Hygienists

Clip-on style Loupes can be attached to virtually any frame of glasses. Flexibility makes this style of loupe very popular with users already wearing glasses, and especially popular with those having prescription lenses

Many users who wear Prescription Glasses are pleased with the Loupes that can be applied to their existing frame, and therefore do not require Custom Made Loupes